## SABLE SLOW SPEED ATTACHMENTS MAINTENANCE AND CLEANING MANUAL

Please retain this manual as your proof of purchase.

MODEL No.	DESCRIPTION	Head Maximum RPM (5K or 20K)
1600201	1:1 NSK Type Contra -Angle	1600401 Snap on Prophyl Head (5K)
1600202	4:1 NSK Type Contra -Angle	1600402 Screw-in Prophyl Head (5K)
1600203	16:1 NSK Type Contra -Angle	1600403 Standard latch Head (5K)
1600204	1:1 MW Type Contra Angle	1600407 BB swing latch Head (20K)
1600210	1:1 NSK Type Contra -Angle	1600405 Miniature latch Head (5K)
1600301	1:1 Straight Nosecone	1600404 BB push Button Head(20K)
1600302	1:1 MW Shorty & Rhino Type Nosecone	1600411 BB BP Latch Head (20K)
1600303	4:1 Straight Nosecone	1600413 Sable friction grip Head (5K)
1600304	1:1 Lab Nosecone	1600412 StarType BB PB latch head (20K)
		1600406 StarType 12 tooth Standard latch head(5K)
1600310	1:1 Straight Nosecone	1600408 StarType snap-on head (5K)
1600501	Midwest type snapon head (5K)	1600409 Star type screw-in head(5K)
1600502	Midwest type screwin head (5K)	1600410 Star type BB SL head (20K)
1600503	Midwest Standard latch head (5K)	1600504 Midwest BB PL head (5K)

### SAFETY PRECAUTIONS:

1. Always exercise extreme caution when operating your Slowspeed Attachment
2. Always inspect your Slowspeed Attachment thoroughly before each use.
3. Ensure proper burs are used for specific product. Use only burs that are DIN EN ISO 1797-1 norm and are made from hardened steel.
4. Do not use burs that are bent, worn or show any signs of damaged.
5. Ensure that the proper speed(rpm) attachment is used when selecting attachment and motor for Heads
6. Always use PPE (Proper Protective Equipment) when operating and cleaning your Slowspeed Attachment See your Occupational Health and Safety Act or other safety ordinances for specific instructions.
7. If product is dropped your dealer should be contacted, it is recommended that you do not use.
8. Contact your dealer right away if any problem should occur with any of your Sable products.
9. Sterilize before first use and before each subsequent use there after
10. This product is for use by a dentist or trained professional in the use of this product.

### FEATURES

1. Attachments conform to International Standard Organization (ISO) specifications.
2. Will attach to all E type (ISO) Air Motors.
3. Autoclavable.

### SPECIFICATIONS

1. Connection with Air Motor: ISO Standard
2. Method of Sterilization: Autoclave

### HOW TO USE

1. Slide Attachment onto shaft of E Type or MW type Air Motor until locked in place.
2. To remove Attachment: pull Attachment, pushing release button located on the exterior of Air Motor.
3. Direction of operation is determined by Air Motor.

### LUBRICATION

USE ONLY AUTOCLAVABLE HANDPIECE LUBRICANT.

1. Remove Attachment from Air Motor.
2. Remove Head from Attachment.
3. Apply lubricant into both ends of Attachment and reapply to Air Motor.
4. Run Air Motor and Attachment at normal speed for 10 seconds.
5. Wipe off excess oil from Attachment surface.

### AUTOCLAVING PROCEDURE

1. Wipe off debris on the Attachment surface with alcohol-soaked cloth.
2. Lubricate the Attachment. (Refer to Lubrication Instructions).
3. Autoclave for 20 minutes at 135 C.

### ATTENTION

This is a precision instrument. Do not drop or apply excessive force when handling.

### WARRANTY INFORMATION

The Sable slow speed Attachment is warranted to be free of manufacturing defects for a period of 3 months\* from date of purchase.

\*Warranty is conditional on proper use and maintenance of Air Motor. Warranty is limited to repair or replacement at manufacturer's discretion. Under no circumstances will the manufacturer be held responsible for incremental or consequential damages.

Please contact point of purchase authorized distributor for warranty repair.

Sable Industries Inc.  
100 Campbell Ave. Suite 5  
Kitchener, Ontario, Canada  
N2H 4X8

File no: PPIS1005 Rev-E,  
SableSlowSpeedEtyp1  
Date: 04/30/2014  
PRINTED IN CANADA

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